
Actim PROM Test (IGFBP-1)

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Evaluation studies on the Actim PROM test as a method to detect premature rupture of fetal membranes

Publication	Number of patients	Gestational age (wk)	End-point	Sensitivity	Specificity	NPV	PPV
Rutanen et al. 1996	130	15-37	Clinical confirmation	100%	94.7%	100%	93.2%
Ragosch et al. 1996	75	22-41	Clinical confirmation, Dye injection	100%	83%	100%	83%
Jain et al. 1998	100	24-42	Clinical confirmation	100%	89%	100%	75%
Erdemoglu et al. 2004	71	31.9±5.3	Clinical confirmation	97%	97%	97%	97%
Kubota et al. 1998	48	15-41	Clinical confirmation	94.7%	93.1%	N/A	N/A

Actim PROM is a highly specific and sensitive method to detect a membrane rupture.

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Comparison of methods to detect PROM

Kubota & Takeuchi 1998

	Actim PROM	pH	Ferning
Sensitivity	94.7	73.3	42.1
Specificity	93.1	72.4	75.9

Actim PROM is the most reliable test for detecting premature rupture of fetal membranes.

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IGFBP-1 concentration in various body fluids

Rutanen et al. 1993

Sample	Concentration of IGFBP-1
Normal adult serum	0.5–30 $\mu\text{g/l}$
Serum (pregnancy)	58–600 $\mu\text{g/l}$
Urine	Undetectable
Semen	Undetectable
Amniotic fluid	10 000–400 000 $\mu\text{g/l}$

IGFBP-1 concentration is extremely high in amniotic fluid only.

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References

Erdemoglu E and Mungan T. Significance of detecting insulin-like growth factor binding protein-1 in cervicovaginal secretions: Comparison with nitrazine test and amniotic fluid volume assessment. Acta Obstet Gynecol Scand (2004) 83: 622-626.

Rutanen E-M. Insulin-like growth factors in obstetrics. Curr Opin Obstet Gynecol (2000) 12: 163-168.

Guibourdenche J et al. Rapid detection of insulin-like growth factor-binding protein-1 and foetal fibronectin in cervico-vaginal secretions to diagnose premature membrane rupture. Ann Clin Biochem (1999) 36: 388-390.

Jain K and Morris P G. A clinical study to evaluate the usefulness of the MAST test in diagnosing pre-labour rupture of membranes. J Obstet Gynaecol (1998) 18: 33-36.

Kubota T and Takeuchi H. Evaluation of insulin-like growth factor binding protein-1 as a diagnostic tool for rupture of the membranes. J Obstet Gynecol Res (1998) 24: 411-417.

Ragosch V et al. Insulin like growth factor binding protein 1 (IGFBP-1) and fetales Fibronectin in der Diagnostik eines vorzeitigen Blasensprunges. GebFra (1996) 56: 1-6.

Rutanen E-M et al. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta (1996) 253: 91-101.

Rutanen E-M et al. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check Membrane immunoassay in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta (1993) 214: 73-81.

Wathen et al. Levels of insulin-like growth factor-binding protein-1 increase rapidly in amniotic fluid from 11 to 16 weeks of pregnancy. J Endocrinol (1993) 137:R1-R4.

Patents

EP0677170, EP0565541, US5554504, US5712170, US5965458

Ordering information

Product Description	REF number
Actim PROM 20 test kit	30832ETAC
Actim PROM 10 test kit	30831ETAC
Actim PROM 3 test sample kit	30833ETAC
Actim PROM Controls	30800ETAC